

Hybridization Protocol

Product Name: **Matrix II Slides**

Catalog No.: CPX 5, CPX 25, CPX 50

This protocol has been optimized to allow users to produce arrays of cDNA or oligonucleotides dissolved in 50% DMSO, 3x SSC, or 150 mM sodium phosphate buffer.

1. Preparing targets for spotting

- A. We recommend dissolving cDNA, unmodified oligos or amino-modified oligos in 50% DMSO/50% H₂O, 3x SSC, or 150 mM sodium phosphate buffer.
- B. Transfer the targets to a 96 or 384 spotting plate with a volume of ~20 μ L/well.
- C. Gently shake the plate to bring liquid to the bottom of the wells.
- D. Now you can set up the array spotter and start printing slides.

2. UV Cross-linking

- A. After spotting, allow the slides to dry for 30 minutes at room temperature.
- B. UV cross-link the slides at **400mJ when working with cDNA**, or **600mJ when working with oligos**.
Important: Do not omit the step!

3. Pre-Treatment

- A. Prepare and preheat Pre-treatment Solution (2x SSC/0.2% SDS/0.1% BSA) to 55 °C.
- B. Place the slides in a slide rack.
- C. Place the slides in a staining jar with Pre-treatment Solution for 20 to 30 minutes at room temperature on an orbital shaker.
- D. Remove the slides from the pre-treatment solution and **rinse thoroughly** with Milli-Q water.
- D. Dry the slides with a gentle stream of nitrogen.
- E. Then allow the slides to dry at room temperature for 15 – 30 minutes.

4. Preparation of probes (Be sure to suspend the probes in FMB Hybridization Buffers)

- A. We recommend using 35 μ L of hybridization mixture for each full slide.
- B. Spin-dry cDNA or oligo probes in a Speedvac.
- C. For hybridization of each slide, resuspend the probes in 3 μ L of nuclease free water.
- D. Quickly vortex the probes and centrifuge it for 30 seconds.
- E. Denature the probes on a heat-block at 90 – 95°C for 3 – 5 minutes.
- F. Remove the probe mix from the heat-block and immediately place it on ice.
- G. Add 32 μ L of hybridization buffer to probes. Please use FMB cDNA Hybridization Buffer (P/N: HBC 01) when working with cDNA or FMB Oligo Hybridization Buffer (P/N: HBO 01) when working with oligos.
- H. Place the mixture on ice for five minutes before applying it to slides.

5. Hybridization

- A. We recommend using 35 μ L of the probe mix for each slide when using a full coverslip (24 mm X 60 mm).
- B. Clean coverslips with 70% ethanol, and blow-dry with nitrogen.
- C. Quickly vortex the probe mixture before applying it to the printed slides.
- D. Place 35 μ L of the probe mixture on each slide and carefully lay down the coverslip (Be careful – avoid air bubbles under the coverslip).
- E. Incubate the slides in a humidified chamber with 100% humidity, at 42 °C or the temperature suitable for your samples for 12 – 14 hours.

6. Washing the slides after hybridization

- A. Preheat Wash Solution 1 (0.2x SSC, 0.2% SDS) and Wash Solution 2 (0.2x SSC) to 55 °C.
- B. Remove cover slips by quickly rinsing the slides with running DI water, and then place the slides in a slide rack.
- C. Immerse the slides in Wash Solution 1 for 20 minutes at room temperature on an orbital shaker.
- D. Transfer the slides to a staining dish with Wash Solution 2 and gently dip the slides up and down for one minute.
- E. Repeat twice. Be sure to use **fresh** Wash Solution 2 each time.
- F. **Thoroughly rinse the slides three times with fresh Milli-Q water** at room temperature.
- G. Dry slides with a gentle stream of nitrogen immediately.
- H. The slides are now ready for scanning.