

## Monoclonal Anti-human Alpha Fetoprotein (AFP)

### Specification Sheet

<b>Catalog #:</b>	MA013
<b>Lot #:</b>	1188
<b>Description:</b>	Monoclonal Anti-Human Alpha Fetoprotein
<b>Specificity:</b>	This antibody reacts with AFP. There is no cross reactivity with human albumin, CEA, ferritin, hCG and PSA.
<b>Clone:</b>	2B51G8
<b>Host Animal:</b>	Hybridization of P3x63-Ag 9.653 myeloma cells with spleen cells from Balb/c mouse
<b>Isotype:</b>	IgG1
<b>Source:</b>	Mouse ascites
<b>Immunogen:</b>	AFP purified from human cold blood
<b>Format:</b>	Purified, liquid
<b>Purification:</b>	DEAE chromatography
<b>Concentration:</b>	7.37 mg/mL (OD280 nm, E <sup>0.1%</sup> = 1.40)
<b>Affinity Constant:</b>	Not determined
<b>Buffer:</b>	0.015M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN <sub>3</sub> pH=7.2
<b>Pairing Assays:</b>	Pairs with MA014
<b>Purity:</b>	>95% by HPLC and SDS-PAGE
<b>Application:</b>	AFP quantitative assays by EIA
<b>Storage:</b>	2-8° C for 1 year. Store at -20° C for long term.

**For research use only**