

Monoclonal Anti-human Prostate Specific Antigen (PSA)

Specification Sheet

Catalog #:	MP016
Lot #:	D3008
Description:	Monoclonal Anti-Human PSA-ACT
Specificity:	This antibody reacts with PSA-ACT complex. There is no cross reactivity with human kallikrein and ACT.
Clone:	5B4
Host Animal:	Hybridization of P3x63-Ag 9.653 myeloma cells with spleen cells from Balb/c mouse
Isotype:	IgG1
Source:	Mouse ascites
Immunogen:	PSA purified from human seminal fluid
Format:	Purified, liquid
Purification:	DEAE chromatography
Concentration:	5.84 mg/mL (OD280 nm, E ^{0.1%} = 1.40)
Affinity Constant:	Not determined
Buffer:	0.015M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN ₃ pH=7.2
Pairing Assays:	Pairs with MP010, MP011, MP013, MP015 and MP019
Purity:	>95% by HPLC and SDS-PAGE
Application:	PSA-ACT quantitative assays by EIA
Storage:	2-8° C for 1 year. Store at -20° C for long term.

For research use only