


**BTK (Phospho-Tyr223) Antibody**

<b>Catalog No.</b>	65675
<b>Product Name</b>	BTK (Phospho-Tyr223) Antibody
<b>Size</b>	50ug, 100ug
<b>Concentration</b>	1 mg/ml
<b>Specificity</b>	Detects endogenous levels of BTK only when phosphorylated at tyrosine 223.
<b>Reactivity</b>	Human, Mouse
<b>Immunogen</b>	Synthesized non-phosphopeptide derived from human BTK around the phosphorylation site of tyrosine 223 (A-L-Y <sup>P</sup> -D-Y).
<b>Applications</b>	WB: 1/500 – 1/1000
<b>Swiss-Prot No.</b>	Q06187 (Human)
<b>Host Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Storage</b>	-20°C for one year. Avoid freeze / thaw cycles.
<b>Purification</b>	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific phosphopeptide. The antibody against non-phosphopeptide was removed by chromatography using non-phosphopeptide corresponding to the phosphorylation site.

**Images****For Research Use Only**