


ATF-2 (Ab-73 or 55) Antibody

Catalog No.	75353
Product Name	ATF-2 (Ab-73 or 55) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of total ATF-2 protein.
Reactivity	Human, Mouse, Rat
Immunogen	Synthesized phosphopeptide derived from human ATF-2 around the phosphorylation site of threonine 73 or 55 (T-P-T ^P -R-F).
Applications	WB: 1/500 – 1/1000; IHC: 1/50 – 1/100
Swiss-Prot No.	P15336 (Human); P16951 (Mouse); Q00969 (Rat)
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Storage	-20°C for one year. Avoid freeze / thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images



For Research Use Only