



STAT5A (Ab-694) Antibody

| | |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No. | 75369 |
| Product Name | STAT5A (Ab-694) Antibody |
| Size | 50ug, 100ug |
| Concentration | 1 mg/ml |
| Specificity | Detects endogenous levels of total STAT5 protein. |
| Reactivity | Human, Mouse, Rat |
| Immunogen | Synthesized non-phosphopeptide derived from human STAT5A around the phosphorylation site of tyrosine 694 (D-G-Y ^P -V-K). |
| Applications | WB: 1/500 – 1/1000; IHC: 1/50 – 1/100 |
| Swiss-Prot No. | P42229 (Human) |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Buffer | Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | -20°C for one year. Avoid freeze/thaw cycles. |
| Purification | Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |

Images



For Research Use Only