



BCL-2(Ab-56) Antibody

Catalog No.	75380
Product Name	BCL-2(Ab-56) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of total BCL-2 protein.
Reactivity	Human
Immunogen	Synthesized non-phosphopeptide derived from human BCL-2 around the phosphorylation site of threonine 56 (G-H-T ^P -P-H).
Tested Applications	WB: 1/500 – 1/1000; ELISA
Entrez Gene No.	596 (Human)
Swiss-Prot No.	P10415 (Human)
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Storage	-20°C for one year. Avoid freeze/thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images



Peptide - +

Immunohistochemical analysis of paraffin- embedded human breast carcinoma tissue using BCL-2(Ab-56) antibody.

For Research Use Only