


BCL-XL (Ab-62) Antibody

Catalog No.	75382
Product Name	BCL-XL (Ab-62) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of total BCL-XL protein.
Reactivity	Human, Mouse, Rat
Immunogen	Synthesized non-phosphopeptide derived from human BCL-XL around the phosphorylation site of serine 62 (A-D-S ^P -P-A).
Applications	Western Blot: 1/500 – 1/1000
Entrez Gene No.	598 (Human); 12048 (Mouse); 24888 (Rat)
Swiss-Prot No.	Q5TE63 (Human); Q07817 (Human); Q64373 (Mouse); Q52KS0 (Rat)
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Storage	-20°C for one year. Avoid freeze / thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images**For Research Use Only**