


Tau(Ab-205) Antibody

Catalog No.	75418
Product Name	Tau(Ab-205) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of total tau protein.
Reactivity	Human, Mouse, Rat
Immunogen	Synthesized non-phosphopeptide derived from human Tau around the phosphorylation site of threonine 205 (P-G-T ^P -P-G).
Applications	WB: 1/500 – 1/1000
Swiss-Prot No.	P10636 (Human)
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Buffer	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 240mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	-20°C for one year. Avoid freeze/thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images



Western blot analysis of extract from mouse brain tissue using Tau (Ab-205) antibody (Lane 1 and 2) and Tau (phospho- Thr205) antibody (Lane 3 and 4).

For Research Use Only