


ASK1 (Ab-83) Antibody

Catalog No. 75446
Product Name ASK1 (Phospho-Ser83) Antibody
Size 50ug, 100ug
Concentration 1 mg/ml
Specificity Detects endogenous levels of total ASK1 protein.
Reactivity Human, Mouse
Immunogen Synthesized phosphopeptide derived from human ASK1 around the phosphorylation site of serine 83 (G-S-S^P-V-G).
Applications WB: 1/500 – 1/1000; IHC: 1/50 – 1/100
Swiss-Prot No. Q98383 (Human)
Host Species Rabbit
Clonality Polyclonal
Isotype IgG
Storage -20°C for one year. Avoid freeze / thaw cycles.
Purification Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images



For Research Use Only