



SAPK/JNK (Ab-183) Antibody

Catalog No.	75562
Product Name	SAPK/JNK (Ab-183) Antibody
Size	50ug, 100ug
Concentration	1 mg/ml
Specificity	Detects endogenous levels of total SAPK/JNK protein.
Reactivity	Human, Mouse, Rat
Immunogen	Synthesized non-phosphopeptide derived from human SAPK/JNK around the phosphorylation site of threonine 183 (M-M-T ^P -P-Y).
Applications	WB: 1/500 – 1/1000; IHC: 1/50 – 1/100
Swiss-Prot No.	P45984; P53779
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Buffer	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	-20°C for one year. Avoid freeze/thaw cycles.
Purification	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Images



Western blot analysis of extracts from 293 cell using SAPK/JNK (Ab-183) Antibody (Lane 1, 2) and SAPK/JNK (phospho-Thr183) antibody (Lane3 and 4).



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using SAPK/JNK (Ab-183) antibody.

Left: peptide -
Right: peptide +

For Research Use Only