


**SAPK/JNK (Ab-185) Antibody**

<b>Catalog No.</b>	75563
<b>Product Name</b>	SAPK/JNK (Ab-185) Antibody
<b>Size</b>	50ug, 100ug
<b>Concentration</b>	1 mg/ml
<b>Specificity</b>	Detects endogenous levels of total SAPK/JNK protein.
<b>Reactivity</b>	Human, Mouse, Rat
<b>Immunogen</b>	Synthesized non-phosphopeptide derived from human SAPK/JNK around the phosphorylation site of tyrosine 185 (T-P-Y <sup>P</sup> -V-V).
<b>Applications</b>	IHC: 1/50 – 1/100
<b>Swiss-Prot No.</b>	P45984; P53779
<b>Host Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Buffer</b>	Phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Storage</b>	-20°C for one year. Avoid freeze/thaw cycles.
<b>Purification</b>	Affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

**Images**

Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using SAPK/JNK (Ab-185) antibody.

Left: peptide –  
Right: peptide +

**For Research Use Only**