

Monoclonal Anti-human Alpha Fetoprotein (AFP)

Specification Sheet

Catalog #:	MA012
Lot #:	1210
Description:	Monoclonal Anti-Human Alpha Fetoprotein
Specificity:	This antibody reacts with AFP. There is no cross reactivity with human albumin, CEA, ferritin, hCG and PSA.
Clone:	3A1
Host Animal:	Hybridization of P3x63-Ag 9.653 myeloma cells with spleen cells from Balb/c mouse
Isotype:	IgG1
Source:	Mouse ascites
Immunogen:	AFP purified from human cold blood
Format:	Purified, liquid
Purification:	DEAE chromatography
Concentration:	5.12 mg/mL (OD280 nm, E ^{0.1%} = 1.40)
Affinity Constant:	Not determined
Buffer:	0.015M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN ₃ pH=7.2
Pairing Assays:	Pairs with MA011
Purity:	>95% by HPLC and SDS-PAGE
Application:	AFP quantitative assays by EIA
Storage:	2-8° C for 1 year. Store at -20° C for long term.

For research use only