

## Monoclonal Anti-Estriol

### Specification Sheet

<b>Catalog #:</b>	ME041
<b>Description:</b>	Monoclonal Antibody to Estriol
<b>Specificity:</b>	Estriol            100% Estradiol        <0.1%
<b>Clone:</b>	1D92D5
<b>Host Animal:</b>	Hybridization of P3x63-Ag 9.653 myeloma cells with spleen cells from Balb/c mouse
<b>Isotype:</b>	IgG1
<b>Source:</b>	Mouse ascites
<b>Immunogen:</b>	Estriol-BSA
<b>Format:</b>	Purified, liquid
<b>Purification:</b>	DEAE chromatography
<b>Affinity Constant:</b>	Not determined
<b>Buffer:</b>	0.015M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN <sub>3</sub> , pH=7.2
<b>Purity:</b>	>95% by HPLC and SDS-PAGE
<b>Concentration:</b>	4.13 mg/ml, Lot D4216 (OD280 nm, E <sup>0.1%</sup> = 1.40)
<b>Pairing Assays:</b>	N/A
<b>Application:</b>	Estriol quantitative assays by EIA
<b>Storage:</b>	2-8° C for 1 year. Store at -20° C for long term.

**For research use only**