

## Monoclonal Anti-human Laminin

### Specification Sheet

<b>Catalog #:</b>	ML021
<b>Lot #:</b>	1232
<b>Description:</b>	Monoclonal Anti-Human Laminin
<b>Specificity:</b>	This antibody is specific to human laminin. It does not reactive with mouse laminin.
<b>Clone:</b>	1A9E3
<b>Host Animal:</b>	Hybridization of P3x63-Ag 9.653 myeloma cells with spleen cells from Balb/c mouse
<b>Isotype:</b>	IgG1
<b>Source:</b>	Mouse ascites
<b>Immunogen:</b>	Human laminin
<b>Format:</b>	Purified, liquid
<b>Purification:</b>	DEAE chromatography
<b>Concentration:</b>	4.42 mg/mL (OD280 nm, $E^{0.1\%} = 1.40$ )
<b>Affinity Constant:</b>	Not determined
<b>Buffer:</b>	0.015M potassium phosphate buffer, 0.85% NaCl, 0.05% $\text{NaN}_3$ , pH=7.2
<b>Pairing Assays:</b>	Pairs with ML022 and ML023
<b>Purity:</b>	>95% by HPLC and SDS-PAGE
<b>Application:</b>	Human laminin quantitative assays by EIA
<b>Storage:</b>	2-8° C for 1 year. Store at -20° C for long term.

**For research use only**